


GUIDE TO PROJECT CODE (Pre-Processing)

1.0 Setup

- Create a folder to contain your data:




Make sure your subjects have the BOLD and T1 subdirectories with all the data (shown for Sub09).

- Start SPM (inside project folder):




2.0 Creating the Pre-Processing Script


- Click on Batch button



- Add a couple “Named File Selectors” for ANAT, BOLDs




- Setup ANAT



Name it “ANAT”

Choose <Subject>/T1/mprage.nii for “File Set”

- Setup BOLDs








Name it "BOLD"

Choose ALL FILES in desired Subject Run <Subject>/BOLD/<Run_#>/<SELECT ALL>

There should be 208 BOLD images

- Add all pre-processing steps





- Figure out Acquisition parameters

Find Data Acquisition Information from article: <https://www.nature.com/articles/sdata20151>

MRI acquisition


The MRI data were collected from a Siemens 3T TIM TRIO (Siemens, Erlangen, Germany). A standard 1 mm isotropic T1-weighted 'structural' image was acquired using an MPAGE sequence (TR 2,250 ms, TE 2.98 ms, TI 900 ms, 190 Hz/pixel; flip angle 9°). The face of the participant in the T1 image was subsequently manually removed to help maintain anonymity. The functional data were acquired using an EPI sequence of 33, 3 mm-thick axial slices (TR 2000 ms, TE 30 ms, flip angle 78°). Slices were acquired in an interleaved fashion, with odd then even numbered slices (where slice 1 was the most inferior slice) and a 25% distance spacing (increased where necessary to cover whole of cortex), resulting in a range of voxel sizes of 3×3×3.75 mm to 3×3×4.05 mm across participants. 210 volumes were acquired in each of 9 runs (note that the 3 initial TRs were discarded to allow saturation of T1 effects).



Slices: 33



TR: 2


Slice Order: Interleaved Bottom-Up (1 to 33 by 2s, 2 to 32 by 2s)

- Setup Slice Timing













- Setup Realign




- Setup Coregister









- Setup Normalize






- Setup Smoothing




3.0 Running Script

You can run from the Play button for 1 Subject 1 run, but you have to change BOLD named file selector for every Run... and then ANAT and BOLD for every Subject...



OR save the script as Matlab script and make some edits:



Edits to Script – **Add init section to start of Script:**

Yellow – Change for your subject/run

Orange – Sets up the ANAT file for the subject

Red – Sets up the BOLD files for subject/run

The **if** statement below BOLD files checks if more files than expected have been imported, and exits script if so. If you want to redo pre-processing, I recommend deleting the Run0# folder and re-copying it from the original data.

```
1 %-----
2 % Job saved on 12-Apr-2025 10:30:05 by cfg_util (rev $Rev: 7345 $)
3 % spm SPM - SPM12 (7771)
4 % cfg_basicio BasicIO - Unknown
5 %-----
6
7 %% INIT
8 % spm('defaults', 'FMRI'); % uncomment if you want to launch script without opening SPM
9 spm_jobman('initcfg');
10 matlabbatch = {};
11
12 % modify your dataset folder here:
13 project_folder = 'C:\Users\samsa\Documents\MATLAB\elec445\project\';
14 % modify your subject here:
15 subject = 'Sub09'; % 1 to 16, keep 0 padding for 1 digit subjects
16 % modify the run here:
17 run = 1; % 1 to 9
18
19 subject_folder = append(project_folder, subject);
20 run_number = append('Run_0', int2str(run));
21
22 ANAT_file = append(subject_folder, '\T1\mprage.nii');
23
24 BOLDS = dir(fullfile(subject_folder, 'BOLD', run_number, '*.nii'));
25 BOLD_files = cell(size(BOLDS, 1), 1);
26 for iFile = 1:size(BOLDS, 1)
27     BOLD_files{iFile} = fullfile(BOLDS(iFile).folder, BOLD(iFile).name);
28 end
29
30 % Check total amount of BOLDS is 208 images
31 if iFile ~= 208
32     disp('There seems to be a problem with the number of BOLD files: ');
33     disp(' - Expected BOLD files: 208');
34     disp([' - Imported BOLD files: ', int2str(iFile)]);
35     disp(['Maybe you already processed this run? ( ', subject, ' : ', run_number, ' )']);
36     return;
37 end
38
39 disp(['>> Pre-Processing: ', subject, ' - ', run_number, '...']);
40
41 %% SETUP FILES
42 matlabbatch{1}.cfg_basicio.file_dir.file_ops.cfg_named_file.name = 'ANAT';
43 matlabbatch{1}.cfg_basicio.file_dir.file_ops.cfg_named_file.files = {{ ANAT_file }};
44 matlabbatch{2}.cfg_basicio.file_dir.file_ops.cfg_named_file.name = 'BOLD';
45 matlabbatch{2}.cfg_basicio.file_dir.file_ops.cfg_named_file.files = { BOLD_files };
46
```

Added call to spm_jobman at the end to run the job:

```
%% RUN PRE-PROCESSING
spm_jobman('serial', matlabbatch)
disp('>> Preprocessing Done!')
```

Final Pre-Processing Script (save as .m file in project folder alongside all Subjects):

```
%-----  
% Job saved on 12-Apr-2025 10:30:05 by cfg_util (rev $Rev: 7345 $)  
% spm SPM - SPM12 (7771)  
% cfg_basicio BasicIO - Unknown  
%-----  
  
%% INIT  
% spm_jobman('serial', matlabbatch) % uncomment if you want to launch script  
without opening SPM  
spm_jobman('initcfg');  
matlabbatch = {};  
  
% modify your dataset folder here:  
project_folder = 'C:\Users\samso\Documents\MATLAB\elec445\project\';  
% modify your subject here:  
subject = 'Sub09'; % 1 to 16, keep 0 padding for 1 digit subjects  
% modify the run here:  
run = 1; % 1 to 9  
  
subject_folder = append(project_folder, subject);  
run_number = append('Run_0', int2str(run));  
  
ANAT_file = append(subject_folder, '\T1\mprage.nii');  
  
BOLDS = dir(fullfile(subject_folder, 'BOLD', run_number, '*.nii'));  
BOLD_files = cell(size(BOLDS, 1), 1);  
for iFile = 1:size(BOLDS, 1)  
    BOLD_files{iFile} = fullfile(BOLDS(iFile).folder, BOLDS(iFile).name);  
end  
  
if iFile ~= 208  
    disp('There seems to be a problem with the number of BOLD files: ');  
    disp(' - Expected BOLD files: 208');  
    disp([' - Imported BOLD files: ', int2str(iFile)]);  
    disp(['Maybe you already processed this run? ( ', subject, ' : ', run_number, '  
)']);  
    return;  
end  
  
disp(['>> Pre-Processing: ', subject, ' - ', run_number, '...']);  
  
%% SETUP FILES  
matlabbatch{1}.cfg_basicio.file_dir.file_ops.cfg_named_file.name = 'ANAT';  
matlabbatch{1}.cfg_basicio.file_dir.file_ops.cfg_named_file.files = {{ ANAT_file  
}};  
matlabbatch{2}.cfg_basicio.file_dir.file_ops.cfg_named_file.name = 'BOLD';  
matlabbatch{2}.cfg_basicio.file_dir.file_ops.cfg_named_file.files = { BOLD_files  
}';  
  
%% SLICE TIMING  
matlabbatch{3}.spm.temporal.st.scans{1}(1) = cfg_dep('Named File Selector: BOLD(1)  
- Files', substruct('.', 'val', '{}', {2}, '.', 'val', '{}', {1}, '.', 'val', '{}', {1},  
'.', 'val', '{}', {1}), substruct('.', 'files', '{}', {1}));
```

```

matlabbatch{3}.spm.temporal.st.nslices = 33;
matlabbatch{3}.spm.temporal.st.tr = 2;
matlabbatch{3}.spm.temporal.st.ta = 1.93939393939394;
matlabbatch{3}.spm.temporal.st.so = [1 3 5 7 9 11 13 15 17 19 21 23 25 27 29 31 33
2 4 6 8 10 12 14 16 18 20 22 24 26 28 30 32];
matlabbatch{3}.spm.temporal.st.refslice = 16;
matlabbatch{3}.spm.temporal.st.prefix = 'a';

```

%% REALIGN

```

matlabbatch{4}.spm.spatial.realign.estwrite.data{1}(1) = cfg_dep('Slice Timing:
Slice Timing Corr. Images (Sess 1)', substruct('.', 'val', '{}', {3}, '.', 'val',
 '{}', {1}, '.', 'val', '{}', {1}), substruct('()', {1}, '.', 'files'));
matlabbatch{4}.spm.spatial.realign.estwrite.eoptions.quality = 0.9;
matlabbatch{4}.spm.spatial.realign.estwrite.eoptions.sep = 4;
matlabbatch{4}.spm.spatial.realign.estwrite.eoptions.fwhm = 5;
matlabbatch{4}.spm.spatial.realign.estwrite.eoptions.rtm = 1;
matlabbatch{4}.spm.spatial.realign.estwrite.eoptions.interp = 2;
matlabbatch{4}.spm.spatial.realign.estwrite.eoptions.wrap = [0 0 0];
matlabbatch{4}.spm.spatial.realign.estwrite.eoptions.weight = '';
matlabbatch{4}.spm.spatial.realign.estwrite.roptions.which = [2 1];
matlabbatch{4}.spm.spatial.realign.estwrite.roptions.interp = 4;
matlabbatch{4}.spm.spatial.realign.estwrite.roptions.wrap = [0 0 0];
matlabbatch{4}.spm.spatial.realign.estwrite.roptions.mask = 1;
matlabbatch{4}.spm.spatial.realign.estwrite.roptions.prefix = 'r';

```

%% COREGISTER

```

matlabbatch{5}.spm.spatial.coreg.estwrite.ref(1) = cfg_dep('Named File Selector:
ANAT(1) - Files', substruct('.', 'val', '{}', {1}, '.', 'val', '{}', {1}, '.', 'val',
 '{}', {1}, '.', 'val', '{}', {1}), substruct('.', 'files', '{}', {1}));
matlabbatch{5}.spm.spatial.coreg.estwrite.source(1) = cfg_dep('Realign: Estimate &
Reslice: Mean Image', substruct('.', 'val', '{}', {4}, '.', 'val', '{}', {1},
 '.', 'val', '{}', {1}, '.', 'val', '{}', {1}), substruct('.', 'rmean'));
matlabbatch{5}.spm.spatial.coreg.estwrite.other(1) = cfg_dep('Realign: Estimate &
Reslice: Resliced Images (Sess 1)', substruct('.', 'val', '{}', {4}, '.', 'val',
 '{}', {1}, '.', 'val', '{}', {1}, '.', 'val', '{}', {1}), substruct('.', 'sess',
 '()', {1}, '.', 'rfiles'));
matlabbatch{5}.spm.spatial.coreg.estwrite.eoptions.cost_fun = 'nmi';
matlabbatch{5}.spm.spatial.coreg.estwrite.eoptions.sep = [4 2];
matlabbatch{5}.spm.spatial.coreg.estwrite.eoptions.tol = [0.02 0.02 0.02 0.001
0.001 0.001 0.01 0.01 0.01 0.001 0.001 0.001];
matlabbatch{5}.spm.spatial.coreg.estwrite.eoptions.fwhm = [7 7];
matlabbatch{5}.spm.spatial.coreg.estwrite.roptions.interp = 4;
matlabbatch{5}.spm.spatial.coreg.estwrite.roptions.wrap = [0 0 0];
matlabbatch{5}.spm.spatial.coreg.estwrite.roptions.mask = 0;
matlabbatch{5}.spm.spatial.coreg.estwrite.roptions.prefix = 'r';

```

%% NORMALISE

```

matlabbatch{6}.spm.spatial.normalise.estwrite.subj.vol(1) = cfg_dep('Named File
Selector: ANAT(1) - Files', substruct('.', 'val', '{}', {1}, '.', 'val', '{}', {1},
 '.', 'val', '{}', {1}, '.', 'val', '{}', {1}), substruct('.', 'files', '{}', {1}));
matlabbatch{6}.spm.spatial.normalise.estwrite.subj.resample(1) =
cfg_dep('Coregister: Estimate & Reslice: Resliced Images', substruct('.', 'val',
 '{}', {5}, '.', 'val', '{}', {1}, '.', 'val', '{}', {1}, '.', 'val', '{}', {1}),
 substruct('.', 'rfiles'));
matlabbatch{6}.spm.spatial.normalise.estwrite.eoptions.biasreg = 0.0001;

```

```


matlabbatch{6}.spm.spatial.normalise.estwrite.eoptions.biasfwhm = 60;
matlabbatch{6}.spm.spatial.normalise.estwrite.eoptions.tpm =
{'C:\Users\samso\Documents\MATLAB\spm12\tpm\TPM.nii'};
matlabbatch{6}.spm.spatial.normalise.estwrite.eoptions.affreg = 'mni';
matlabbatch{6}.spm.spatial.normalise.estwrite.eoptions.reg = [0 0.001 0.5 0.05
0.2];
matlabbatch{6}.spm.spatial.normalise.estwrite.eoptions.fwhm = 0;
matlabbatch{6}.spm.spatial.normalise.estwrite.eoptions.samp = 3;
matlabbatch{6}.spm.spatial.normalise.estwrite.woptions.bb = [-78 -112 -70
78 76 85];
matlabbatch{6}.spm.spatial.normalise.estwrite.woptions.vox = [2 2 2];
matlabbatch{6}.spm.spatial.normalise.estwrite.woptions.interp = 4;
matlabbatch{6}.spm.spatial.normalise.estwrite.woptions.prefix = 'w';

%% SMOOTHING
matlabbatch{7}.spm.spatial.smooth.data(1) = cfg_dep('Normalise: Estimate & Write:
Normalised Images (Subj 1)', substruct('.', 'val', '{}', {6}, '.', 'val', '{}', {1},
.', 'val', '{}', {1}, '.', 'val', '{}', {1}), substruct('()', {1}, '.', 'files'));
matlabbatch{7}.spm.spatial.smooth.fwhm = [8 8 8];
matlabbatch{7}.spm.spatial.smooth.dtype = 0;
matlabbatch{7}.spm.spatial.smooth.im = 0;
matlabbatch{7}.spm.spatial.smooth.prefix = 's';

%% RUN PRE-PROCESSING
spm_jobman('serial', matlabbatch)
disp('>> Preprocessing Done!')

```



Ok, so it takes around ~15 minutes for my CPU to do the pre-processing for 1 Run for a subject.
Make sure the above script is saved in project folder beside Sub folders:



Edit these in script for your Subject/Run:


```
11  
12 % modify your dataset folder here:  
13 project_folder = 'C:\Users\samso\Documents\MATLAB\elec445\project\';  
14 % modify your subject here:  
15 subject = 'Sub09'; % 1 to 16, keep 0 padding for 1 digit subjects  
16 % modify the run here:  
17 run = 1; % 1 to 9  
18
```

Now run using Play button (or type script name in command window):




4.0 Visualizing the Pre-Processing Results


Use the display button to check Slice Timing/Realign/Coregister with T1:




Select the T1 anat as image:




Add an overlay to the T1:




Type 'rra' in filter and select first Coregistered image:




Now make sure it 'overlaps' well:




You can also use the “Check Reg” button (maybe is clearer):




First select the ANAT (**BUT DON'T CLICK DONE YET**):



Now select a file to compare (like first Coregisterer image with 'rra' filter):




Click Done and make sure they match (cursor position on both images):






ALSO MAKE SURE TO CHECK NORMALIZE/SMOOTHED IMAGES

Find Template T1 in your SPM installation folders:



Select a Normalize/Smoothed image (filter 'swrra' in BOLD run for Smoothed):





Use Check Reg or Display, whatever you prefer.